



LEHIGH VALLEY DENTAL REPAIR

SERVING THE DENTAL COMMUNITY FOR OVER 30 YEARS

610-837-8447

15 POINT PREVENTIVE MAINTENANCE PLAN

TO MINIMIZE DOWNTIME

- Check the air compressor oil, filters & moisture.
- Check all unit filters (change if necessary).
- Check all lights for drift (adjust if needed).
- Check x-ray units for drift & voltage and secure mounting to wall.
- Visual inspection of chairs (loose components & fluid leaks).
- Check master water filter (change if needed).
- Check water quality in operatories.
- Inspect vacuum pump for proper function.
- Check amalgam separator filter (replace if needed).
- Check air/water syringe buttons and o rings.
- Check handpiece pressure and tubing's.
- Check foot controls, tighten (lube if needed).
- Check autoclave/statim - gasket and filters.
- Check dental cabinetry for proper function.
- Check Dr. and Assistant stools (wheels and pistons) for proper function.

STARTING AT \$170